



# Healthcare bed accessories

2080 / 3080 / 5384

# Völker accessories


## Healthcare beds 2080, 3080 and 5384.

All Völker accessories can be attached by simply inserting or hooking them into the integrated mounting sleeves, brackets and rails.

The following pages provide a detailed overview of the accessories available for each healthcare bed.

People requiring care have a particular right to respect and regard for their dignity. Caregivers also have a right to equipment that makes their daily work as easy as possible. We think, develop and manufacture with this in mind.

Hand control holder for the 3080.















Smart holders for smart devices  
Our tablet and mobile phone holders  
are particularly convenient accessories  
for residents who use digital devices.



Configuration		Item no.	2080	3080	5384
	<b>Trapeze bar with handle</b> Material: powder coated	PZP-3060	✓	✓	✓
	<b>Trapeze bar with handle</b> Adjustable with one hand Material: powder coated	PZP-3060/2	✓	✓	✓
	<b>Assist rail spacer</b> 1 piece Material: Stainless steel and high-pressure laminate (HPL)	PZK-969/2	-	✓	✓*
	<b>Assist rail spacer cover</b> 1 piece Colour: Brown	PZP-3059	-	✓	✓*
	<b>Assist rail spacer holder</b>	PZP-3070	-	✓	✓*
	<b>Universal belt holder</b> 1 pair				
	Relocatable without tools	PZK-972/3	✓	✓	-
	Fixed in place by screws	PZK-972/4	✓	✓	-
	<b>Crutch holder</b> 1 pair	PZP-965	✓	✓	✓
	<b>Hand control holder</b>	PZK-1044/3	-	✓	-
	<b>Care Clip cable holder</b> For cable routing on the trapeze bar e.g. for a nurse call button	KT3859	✓	✓	✓
	<b>Infusion holder</b> For trapeze bar PZP-3060	PZP-3055/2	✓	✓	✓
	<b>Infusion stand</b> Adjustable with one hand 2 kg load capacity per hook Material: stainless steel *PZP-5084 adapter required	PZP-3057	✓	✓	✓
	<b>ISO standard rail</b> Material: stainless steel				
	Length: 40 cm	PZK-1083/6	✓	✓	-
	Length: 60 cm	PZK-1083/8	✓	✓	-

With some accessories, availability also depends on the year of manufacture and configuration of your bed. Need to know exactly which accessory versions are suitable for your Völkner bed? Then please check the serial number of your bed and contact our customer service.

Configuration		Item no.	2080	3080	5384
	<b>Side rail protective cover</b> For bed length: 200 cm Height of the side rail: 34 cm Colour: Brown	PZP-2056-2.0	✓	-	-
	<b>Side rail protective covers</b> For bed length: 200 cm Height of the side rail: 34 cm Colour: Brown	PZP-3056_C_2.0	-	✓	✓*
	<b>Universal hook</b> 1 pair	PZK-1043	✓	✓	-
	<b>Universal hook</b> 1 pair For ISO standard rails	PZK-943	✓	✓	✓
	<b>Urine bottle basket</b> For hanging on the accessories rail Material: Chrome	PZK-1040	✓	✓	-
	<b>Urine bottle basket</b> Inserts into the mounting sleeve Material: Chrome	PZK-1040/2	-	✓	✓*
	<b>Wall protector wheels</b> Vertical	PZP-3068	✓	✓	✓
	<b>Wall spacer</b> Head side	PMAYA02	-	-	✓
	<b>Wall spacer</b> Long side	PZK-1092	-	-	✓
	<b>Maya lamp</b> Connected behind the headboard using an adapter	PMAYA02	✓	✓	✓
	<b>Tablet holder</b> Can be plugged into the socket of the hand control	PZK-1092	-	✓	✓
	<b>Mobile phone holder</b> Can be plugged into the socket of the hand control	PZK-1093	-	✓	✓

✓ available - not available \* not with fixed side rails

24/7

**ALWAYS  
AVAILABLE**



**FOCUS  
ON SERVICE**



**HIGH-QUALITY  
MATERIALS**



**CERTIFIED  
SUSTAINABILITY**

**Head office**

**Völker GmbH**

Wullener Feld 79  
58454 Witten  
Germany  
+49 (0) 2302 960 960

**Völker Schweiz GmbH**

Moosstrasse 2  
6003 Luzern  
Switzerland  
+41 (0) 44 563 1939

**Völker Medizinprodukte GmbH**

Kaufmannstraße 14  
6020 Innsbruck  
Austria  
+43 (0) 1 205 8365

**Voelker France SAS**

27 Rue Maurice Flandin  
69444 Lyon Cedex 03  
France  
+33 (0) 18 288 1530

**Voelker Nederland B.V.**

Laan 21  
1741 EA Schagen  
Netherlands  
+31 (0) 20 808 1603

**General enquiries**

info@voelker.de  
+49 (0) 2302 960 960

**Sales / Service / Replacement parts**

kundenservice@voelker.de  
+49 (0) 2302 960 9630